

Coming Soon - 2020 1Q

BIO TECHNOLOGY
boditech

ichroma™ ToRC test Panel

Antibody test for Toxoplasma, Rubella, and Cytomegalovirus

- The world first Point-of Care (POC) test for ToRC multiplexing assay
- High sensitivity immuno-fluorescence test
- Three in one combo test: ToRC IgG and ToRC IgM
- Fast output with high accuracy
- For obstetrics and Neonatology
- Versatile platform: ichroma™ II , ichroma™ III and ichroma™ 10






Find out more on

www.boditech.co.kr



Specifications

Items	ToRC IgG	ToRC IgM	Toxoplasma IgG/IgM
Description	Qualitative test of human IgG against Toxoplasma gondii, Cytomegalovirus (CMV), Rubella virus	Qualitative test of human IgM against Toxoplasma gondii, Cytomegalovirus (CMV), Rubella virus	Qualitative test of human IgM and IgG against Toxoplasma gondii
Principle	Fluorescence immunoassay(FIA)		
Sample type	Serum/Plasma/Whole blood		
Reaction time	12 min		
Storage temp	4-30 °C		
Reaction Temp	17-35 °C		
Shelf life	20 months		
Result type	Qualitative but reference value-provided, IU/mL	Qualitative but reference value-provided, COI	Qualitative but reference value-provided, IgM: COI IgG: UI/mL

Items	ToRC IgG	ToRC IgM	Toxoplasma IgG/IgM
Platform			
ichroma™ II	●	●	●
ichroma™ III *	●	●	●
ichroma™ 10 *	●	●	●

* Coming Soon



43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gangwon-do, 24398, Korea
Tel: +82-33-243-1400, Fax: +82-33-243-9373
www.boditech.co.kr sales@boditech.co.kr